

Humigel™ Replacement Kit

- Oasis® guitar humidifiers are designed to hold .6 cc's of Humigel™ crystals. If a portion of the existing crystals are lost due to spillage, it is important to remove all remaining crystals from the humidifier before refilling with fresh crystals. This way, there are always .6 cc's of crystals in the humidifier. Too few crystals will leave liquid water sloshing around inside the humidifier. Too many crystals will over-absorb water reducing the evaporation rate.
- Before re-filling humidifier with replacement crystals, make certain that all previously existing crystals have been washed out.
- Once clean, take the white refill spoon and scoop 4 level measures (.15 cc's each) into the opening of the humidifier. The package of crystals has sufficient volume to allow for 5-6 refills.
- Fill the syringe with distilled water and fill the humidifier slowly. The humidifier holds about 31 cc's, so it will take slightly more than three syringes of water to fill it the first time.
- Humigel™ Crystals take about 5-10 minutes to absorb the water and turn into a gel complex.
- Once completed, you should only need to use one filling of the syringe (10 cc's) to top off the humidifier once a week when you refill it.